

**KITAZATO**<sup>®</sup>



*Pack of 10 pieces. Available in 5 different colours, Red, White, Yellow, Green and Blue*

#### **INTENDED USE**

Cryotop is used for storage of vitrified human oocytes and embryos.

#### **STORAGE INSTRUCTIONS AND STABILITY**

Store in original sterile pack at 15-30 degree Celsius.

Cryotop is stable until the expiration date shown on package label when stored as directed.

Each Cryotop is intended for single use.

Gytech Pty Ltd  
ABN 12 076 599 570



PO Box 76 Armadale North  
Victoria 3143 Australia  
P +61 3 9822 5911  
F +61 3 9822 4911  
E [sales@gytech.com.au](mailto:sales@gytech.com.au)