

## PIPELLE DE CORNIER MARK II



Pipelle de Cornier®Mark II

### Pipelle de Cornier®Mark II



The Pipelle de Cornier® Mark II is a new device based on the original Pipelle de Cornier®. The endometrial cells (cytology specimen) are collected on the rough surface on 2 cm of its tip, and the endometrial tissue (histology specimen) is collected by the plunger's aspiration.

- Systematic screening for endometrial cancer and hyperplasia
- Detection of luteal phase insufficiency
- Monitoring of hormonal treatment
- Menometrorrhagia with or without HRT
- Late menopause
- Endometrial hypertrophy detected by ultrasonography
- Polyps
- Tamoxifen treatment
- Bacteriological culture

Gytech Pty Ltd  
ABN 12 076 599 570



PO Box 76 Armadale North  
Victoria 3143 Australia  
P +61 3 9822 5911  
F +61 3 9822 4911  
E sales@gytech.com.au