

Thomas Medical K-UMI UTERINE MANIPULATOR



- Sterile. Disposable. Ready to use.
- Pinch-Grip Handle allows for quick and effective repositioning during procedures
- Unique pear shaped balloon conforms to anatomy and allows for a secure seal of the cervical opening to prevent injection leaks.
- Only manipulator with a short, smooth, tapered distal tip, reducing the chance for uterine perforations and providing increased patient comfort during placement
- Only manipulator with a flexible (strain-relief), proximal connector to prevent breaks
- This product does not contain latex.

Product Code: TMI 6003

Description: K UMI Uterine Manipulator (Box of 12)

Gytech Pty Ltd
ABN 12 076 599 570

PO BOX 8631
Armadale, Victoria 3143
Australia

P +61 3 9822 5911
E sales@gytech.com.au

